

Pocket-microscope utilisation

S52.80.91A

To use with the standard products

04 Mar 2004

Directions for use

The POCKET MICROSCOPE is a small size, high performance and durable microscope that could be easily carried in the pocket for use at any time.

The outstanding features are:

- very sharp resolving power throughout its field of view ,
- unbelievably bright and light.

To focus the POCKET MICROSCOPE:

1. on the **rectangular** versions:

- hold it so that the notch at the top of the acrylic cylinder will face towards you,
- place the centre of the object lens right above the Crystal
- place your eye on the eye-piece lens,
- then place the Crystal perpendicular to the axe of the pocket microscope and to observe at an angle which can give you the sharpest image.

